

English	3
български	4
Bosanski	5
Čeština	6
Dansk	7
Deutsch	8
Eestlane	9
Español	10
Français	11
Gaeilge	12
Hrvatski	13
Italiano	14
Latvietis	15
Lietuvis	16
Magyar	17
Nederlands	18
Polskie	19
Português	20
Română	21
Slovenščina	22
Slovenský	23
Suomi	24
Svenska	25
Türk	26
Ελληνικά	27
русский	28
Український	29
한국어	30
日本語	31
简体中文	32

REF SK2
REF SK2-CA

Instructions

Intended Use

This device is intended for medical purposes to illuminate body surfaces. It is used for non-invasive visual examination of intact skin and assists clinicians in the assessment and differentiation of skin abnormalities by enhancing the visualization of surface and subsurface structures.

This battery-operated product is designed for external examination only within professional healthcare facilities by medical professionals.

Check the correct operation of the device before use!

Do not use the device if there are visible signs of damage, malfunction, or the device operates in a manner not consistent with what is described in these instructions for use.

In case of a serious incident with the use of this device, notify DermLite immediately and, if required by local regulations, your national health authority.

CAUTION: Do not look directly into the LED light. Patients must close their eyes during examinations.

NOTICE: This Risk Group 1 product emits 365nm and 405nm UV light. Do not look directly at light. The worst case EHV value of this device is 21.66 at a distance of the closest point of human access (skin) and exposure duration of 10s.

CAUTION: Do not use the device in fire or explosive risk area (e.g., oxygen-rich environment).

CAUTION: Do not expose charging terminals to any foreign/extraneous conductive materials.

CAUTION: During extended operation, the patient-contact surface may reach a temperature of up to 42°C (108°F). Avoid prolonged continuous contact with the same area of skin. Use additional caution with patients who may have reduced sensitivity to heat.

CAUTION: This product contains magnets. Do not use in proximity of pacemakers or other magnetically sensitive devices or items.

WARNING: This product can expose you to chemicals including methylene chloride and hexavalent chromium, which are known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.

Electromagnetic Compatibility

This device complies with the EMC Emissions and Immunity level requirements of the standard IEC 60601-1-2:2014. The emission characteristics of this equipment make it suitable for use in professional healthcare environment as well as residential environment (CISPR 11 Class B). This equipment offers adequate protection to radio communication service. In the rare event of interference to the radio communication service, the user might need to take mitigation measures, such as relocating or re-orienting equipment.

WARNING: Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this

equipment and the other equipment should be observed to verify that they are operating normally.

WARNING: Use of accessories other than those provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

WARNING: Portable RF transmitters should be used no closer than 30 cm (12 inches) to any part of the device. Otherwise, degradation of the performance of this equipment could result.

Do not use this device during magnetic resonance imaging (MRI) or in an MRI environment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Any changes/modifications to this equipment not approved by DermLite could void the user's authority to operate the equipment.

IMPORTANT: Before initial use, please charge the internal lithium ion battery. See "Battery and Charging" for details.

Using Your DermLite skeen

The DermLite skeen is a handheld digital dermatoscope with an integrated display, camera, and lighting system for clinical and dermoscopic skin and hair imaging.

Press the Power button (P) to turn it on. Swipe up to unlock and follow any on-screen prompts.

For clinical imaging, hold the device with the display facing you and the camera (C) aimed at the target area.

For dermoscopy, attach the dermoscopic lens (DL) with spacer (S) and gently press the snap-on contact plate (CP) flat against the skin. For non-contact dermoscopy, remove the contact plate.

Use the two action buttons (AB) to switch between polarized, non-polarized, and UV lighting modes, or to capture images. For more features, refer to the FotoFinder software instructions.

To enter standby mode, briefly press the Power button or place the device on the Charging Base (CB). To fully power off, press and hold the Power button for 2-3 seconds and follow the on-screen instructions.

IceCap® Disposable Caps

To support your infection control measures, snap one of the included IceCap (IC) disposable covers over the installed faceplate (CP). To remove and dispose, simply pull it off.

Battery and Charging

After extensive use or time, battery performance may degrade, causing reduced capacity or swelling. Watch for signs such as shortened battery life, unexpected shutdowns, erratic behavior, or sudden power loss. If observed, stop using the device.

To charge, connect the supplied USB cable to the included power adapter and plug it into

either the Charging Port (PP) at the bottom of the Charging Base (CB) or directly into skeen's USB port.

When using the Charging Base, ensure its charging pins make proper contact. During charging, the indicator light (LED) pulses white and turns solid white when fully charged. It may briefly pulse white when checking charge status. If no power is supplied, the LED remains off.

Dermoscopy Lens (DL)

To attach the lens, align its indicator line (IL) with the corresponding mark on the skeen body. If not properly aligned, its magnets will not engage. Simply pull off to detach.

For non-contact imaging, remove the snap-on contact plate (CP) by gently pulling it off. Reattach by snapping it back into place. To remove the Spacer (S), turn it counter-clockwise and pull it off; reattach with a clockwise turn until secure.

The IceCap storage receptacle built into the charging base may also be used to store the dermoscopy lens.

Troubleshooting

Please check www.derm-lite.com for the most current troubleshooting information. If your device requires servicing, visit www.derm-lite.com/service or contact your local DermLite dealer.

Care and Maintenance

CAUTION: No modification of this equipment is allowed.

Your device is designed for trouble-free operation. Repairs shall be made only by qualified service personnel. It is recommended to check for correct operation of the device every six months.

Cleaning

The exterior of your device, except the optical parts, may be wiped clean with isopropyl alcohol (70% vol.) prior to use on a patient. The lens should be treated as high-quality photographic equipment and should be cleaned with standard lens cleaning equipment and protected from harmful chemicals. Do not use abrasive material on any part of the equipment or immerse the device in liquid. Do not autoclave.

Warranty: 2 years for parts and labor.

Disposal

This device contains electronics and a lithium-ion battery that must be separated for disposal and may not be disposed of with general household waste. Please observe local disposal regulations.

This set includes: DermLite skeen (integrated touchscreen, camera, Wi-Fi, USB-C port, charging contacts, and three-button housing), dermoscopy lens (cross-polarized, non-polarized, parallel-polarized and UV illumination, with removable spacer and snap-on faceplate), charging base, USB-C cable, power adapter, 5 IceCap covers, and microfiber cloth.

Technical Description

Visit www.derm-lite.com/technical or contact your local DermLite dealer.

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
Explanation of Symbols

Do not dispose as unsorted municipal waste

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Explanation of Symbols

Product complies with European Medical Device Regulation (2017/745)

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Explanation of Symbols

European authorized representative

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Caution: UV Light. Do not look directly at light.

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Explanation of Symbols

Caution: Federal law (USA) restricts this device to sale by, or on the order of, a licensed healthcare practitioner

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Indicates the authorised representative in Switzerland

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Responsible Person Limited
Summit House, 4-5 Mitchell Street
Edinburgh, EH6 7BD Scotland

E C Rep GmbH Max-
Högger-Strasse 6
8048 Zürich
info@ch-rep.info

E. C. REP Ltd
5 Fitzwilliam Square East,
Dublin 2, D02 R744, Ireland
Tel: +353 1 2 544 944
Email: info@ecrep.ie

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DermLite LLC